



FDR Visionary Suite Specifications

■ X-ray Generator

- Rated output : 50 kW / 65 kW / 80 kW
- Tube voltage : 40 to 150 kV
- Tube current : 10 to 630 mA (50 kW model)
10 to 800 mA (65 kW model)
10 to 1000 mA (80 kW model)
- AEC : Xe detector-type phototimer receiver
combination up to three receivers

■ X-ray Tube Support

- Ceiling fixture : Fixed rail of 4 / 5.5 m
Moving rail of 2 / 2.6 / 3.3 m
- Movement range : Longitudinal 2.95 m (4 m fixed rail)
Longitudinal 4.45 m (5.5 m fixed rail)
Transversal 1.4 m (2 m moving rail)
Transversal 2.0 m (2.6 m moving rail)
Transversal 2.7 m (3.3 m moving rail)
Vertical 1.6 m
- Rotation : Vertical axis $\pm 180^\circ$
Horizontal axis -180° to $+120^\circ$

■ X-ray Tube Unit

- Maximum anode heat content : 400 KHU
- Maximum anode heat dissipation rate : 2200 HU/s
- Focal spot : 0.6 / 1.2 mm

■ Collimator

- Filtration : Inherent filtration 1.1 mmAl eq.
Added filter of Cu 0.1 / 0.2 / 0.3 mm
- Standard accessories : Auto-filter
Line marker
Detent (fitted at the home position)
- Area dosimeter adapter (Option) :
An adapter for dosimeter manufactured
by VACUTEC/PTW

■ Table

- Tabletop size : 810 $\times$ 2350 mm
- Table height : 535 to 930 mm
- Longitudinal range : ± 375 mm
- Transversal range : ± 125 mm
- Bucky moving range : 800 mm
- Max. load : 295 kg
- Standard accessories : Tracking device
Bucky tracking driver

- Options : Compression belt
Side cassette holder
Grip switch
CFRP tabletop
Hand grip
Drip hanger
Rear foot switch

■ Stand

- Distance between Bucky top edge and floor surface :
Manual : 643 to 2143 mm
Motorized : 671 to 2113 mm
- Tilting angle (Function for BR-120T) : -20° to 90°
- Standard accessories : Tracking device
Stop switch
Foot switch
- Options : Hand grip (mounted on top edge of the Bucky)
Hand grip (mounted on back side of the Bucky)
Cassette holder
Front handle
Both side operation
Compression belt
Patient stand (for long view radiography)
Wall mounting option (for BR-120)



FDR D-EVO Advanced C43A

FDR D-EVO Advanced C43A Specifications

- Scintillator : CsI
- Detector external size :
464.5 $\pm$ 1(W) $\times$ 516.7 $\pm$ 1(D) $\times$ 18 $\pm$ 1(H) mm
*excluding convex part of the cable
- Weight : Approx. 4.5 kg (including battery)
- Pixel size : 150 μ m
- Maximum detecting area : 2816 $\times$ 2816 pixels
- Image preview : less than 2 sec
- Cycle time : less than 8 sec

Specifications are subject to change without notice.
All brand names or trademarks are the property of their respective owners.
All products require the regulatory approval of the importing country.
For details on their availability, contact our local representative.
Actual X-ray images are varied by conditions of X-ray system or subjects or other factors.

FUJIFILM
Value from Innovation

FDR Visionary Suite



Tomosynthesis

Energy Subtraction

Long View Imaging

FUJIFILM

FUJIFILM Corporation

26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN
<http://www.fujifilm.com/products/medical/>

Ref. No. XB-1029ER (SK-16-05-F1079-F9711) Printed in Japan ©2015 FUJIFILM Corporation



Freedom and Flexibility in Imaging

Experience a wide range of applications targeted to improve diagnostic capabilities combined with a precise design that facilitates imaging.

With unique image processing functions that enhance quality while reducing the operational dose rate, the FDR Visionary Suite is the next generation in functional X-ray systems — offering ease of operation with minimal impact for patients.

Compatibility with a Broad Range of Cassette DR Panels Ensures Maximum Flexibility

Wide Array of Applications Support Diagnosis

Tomosynthesis

Energy Subtraction


Long View Imaging

Stress-free, Optimized Imaging Workflow

Using the latest technological developments to allow further quality improvements with a low operational dose rate.

1. Using Fujifilm's "ISS method" reading technology to achieve sharper images


In contrast to conventional FPD offerings this system features an indirect conversion FPD using the "ISS method," where the TFT sensor is placed in front of the scintillation layer instead of behind it. With this proprietary method the scattering/dissipation of X-ray signals is significantly reduced, achieving sharper images with a lower X-ray dose.



2. Improved sensitivity in low-density areas using noise reduction circuitry*


Image noise is reduced by using a proprietary noise reduction circuitry. Granularity in low-density areas is improved, helping to boost image quality.

* FDR D-EVO Advanced C43A




3. Image processing technologies producing optimized images

The system features imaging processing technologies that enable desirable image display. These technologies include "Dynamic Visualization," which optimizes images for diagnosis on-screen, and "Flexible Noise Control (FNC) Processing," which reduces granularity by automatically extracting and excluding image noise components.



FNC processing


Multiple Panel Combinations and Variations



It is possible to select the ideal imaging method for each examination and site requirement by selecting from a broad range of panel sizes and types.


[Full Function Model]

With cutting-edge optional technologies such as Tomosynthesis, additional imaging information can be provided, further improving diagnostic capabilities.



[Standard Model]

This model is capable of a wide-range of general radiography imaging, including long view imaging and an Autopositioning option, enabling optimized work flow.



Panel Name/Panel Type	FDR D-EVO Advanced C43A	FDR D-EVO II C43	FDR D-EVO II C35	FDR D-EVO II G43	FDR D-EVO II G35	FDR D-EVO II C24
Scintillator	CsI	CsI	CsI	GOS	GOS	CsI
Size	43×43 cm	43×43 cm	35×43 cm	43×43 cm	35×43 cm	24×30 cm
Applications*	Tomosynthesis	●	-	-	-	-
	Energy Subtraction	●	●	●	●	-
	Long View Imaging	●	●	●	●	-
Cassette Tray	●	●	●	●	●	(Free exposure position type)

* These applications are optional.

Wide Range of Applications That Contribute to Improving Diagnostic Capabilities


Tomosynthesis*

Freedom to reconstruct and display image slices

With this technology the X-ray tube moves through an arc, acquiring a series of images in a single sweep, which can then be reconstructed to create cross-sectional image slices.

Automatic X-ray dose control and background reconstruction

Using the imaging conditions for a single preliminary image as reference, the conditions for Tomosynthesis imaging are set automatically.




High-precision, high-quality imaging to 150µm

By controlling such items as metal artifacts, high-precision imaging down to 150µm is possible.

Energy Subtraction*


Separates images of soft tissue and bone for improved viewing

This technology takes two X-rays, utilizing the difference in X-ray energy absorption to create specific images of soft tissue and bone, etc. The dose rate is changed automatically between shots.



Controlling motion artifacts

Motion artifacts that may occur between exposures are suppressed by multiple resolution alignment processing, allowing for clear images of soft tissue and bone.




Long View Imaging*


Display full-length images of spine or lower limb

This technology uses multiple images taken in one sweep and automatically stitches them to create images of wide area of up to 160cm upright and 120cm recumbent. Disjointed images caused by patient body movements can be automatically corrected through image alignment.*

Upright



Recumbent



Up to 120cm when in recumbent position

*Depending on the degree of disjoint between images it may not be possible to implement automatic image correction.

* These applications are optional.

Smoothly and Surely Ensuring a Stress-free Imaging Environment



① Preparation



② Patient Guidance and Positioning



③ Taking Images



● Completion of room preparation without touching the system
The system features an auto-positioning function that moves the X-ray tube into position automatically. It is possible to pre-set and restore positions from the image guidance menu.

● X-ray stand
A movable scope of 40 to 190cm from the center of the exposure makes it possible to take images of the entire lower limb from the cervical vertebrae down. The exposure platform can be adjusted from -20 to 90 degrees*, making it possible to take images of the head and upper limbs.




*Option

● X-ray table
Using the foot switch and grip switch* it is possible to adjust the height quickly and easily between 53 and 85cm.




*Option



● Easily define the imaging position for each individual patient
With the auto-tracking function the panel and X-ray tube are automatically kept in alignment, making it possible to focus on patient positioning and care. By switching between automatic and manual functions positioning can be simplified, allowing the operator to maintain full control.

● Radiation field linking function
The pre-selected radiation field size for the area to be imaged is automatically set and alignment of the field to the upper or lower portion of the detector is also automatically performed.

● Setting made easy with an LCD touch panel
The touch panel presents image-related information clearly and also enables settings to be changed easily. It is also possible to change the angle of the square LCD panel by 90 degrees to match the direction of the X-ray tube, making it easy to see at all times.



● Change conditions in the X-ray room using the touch panel
All conditions can be changed using the LCD touch panel on the X-ray tube supporting arm, making it possible to set conditions in the X-ray room alone. The changed conditions are relayed in real time to the controller outside the X-ray room.



● “Sound and light” notifies those away from the machine when an X-ray is being taken
“Ready up” and “X-ray in progress” notifications can be clearly understood by sounds and lights on the frame and hand switch. There is a choice of seven colors for the notification lights.

Provision of Easy-to-use Advanced Applications

Work flow ① Tomosynthesis

1 Imaging
Determine the image position and take an X-ray.

2 Collection of images
Imaging conditions are automatically calculated from the pre-shot and up to 60 images are collected.

3 Images are reconstructed into slices
Metal artifacts are suppressed to create a high-precision image.

Work flow ② Image Stitching

1 Imaging
The parameters to obtain an image of wide area are set and an exposure is performed.

2 Collection of images
Multiple images are taken automatically within the pre-set parameters.

3 Automatic stitching
Multiple images are stitched automatically. Disjoints in stitched images caused by patient body movements are also automatically corrected.